

Anti-gp130/CD130 antibody



Catalog Number: 104480

Product name

Anti-gp130/CD130 antibody

Immunogen

[Rat gp130/CD130 \(His Tag\) recombinant protein](#)

Specificity

Rat gp130 / CD130 / IL6ST

Antibody description

Rabbit monoclonal to gp130/CD130

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat gp130 / CD130 / IL6ST (rR gp130 / CD130 / IL6ST; NP_001008725.2; Met4-Glu618) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



gp130 / CD130 / IL6ST Antibody (PE), Rabbit MAb, Flow cytometric analysis

Flow cytometric analysis of Rat IL6ST(CD130) expression on rat splenocytes. SD splenocytes were stained with PE-conjugated anti-Rat IL6ST(CD130). The histogram were derived from gated events with the forward and side light-scatter characteristics of intact cells.