

# SNX3 antibody

Catalog Number: 115480

## Product name

SNX3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SNX3 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, A549 cells, HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 16-20 kDa

## Preparation

This antibody was obtained by immunization of SNX3 recombinant protein (Accession Number: NM\_003795). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:115480(SNX3 antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells using Catalog No:115480(SNX3 Antibody) at dilution of 1:25 and Alexa Fluor 594-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)