

PPARD antibody



Catalog Number: 107480

Product name

PPARD antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PPARD Mouse Monoclonal antibody. Positive WB detected in Colo320 cells, human heart tissue. Observed molecular weight by Western-blot: 54kd

Preparation

This antibody was obtained by immunization of PPARD recombinant protein (Accession Number: NM_177435). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:107480(PPARD antibody) at dilution of 1:1000