

# beta Arrestin 1 antibody



Catalog Number: 108480

## Product name

beta Arrestin 1 antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

## Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

## Antibody description

beta Arrestin 1 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive IF detected in HepG2 cells. Positive IP detected in mouse lung tissue. Positive WB detected in mouse lung tissue, A549 cells, U-937 cells. Observed molecular weight by Western-blot: 55 kDa

## Validations




## Preparation

This antibody was obtained by immunization of beta Arrestin 1 recombinant protein (Accession Number: NM\_004041). Purification method: Antigen affinity purified.

mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:117125(ARRB1 antibody) at dilution of 1:500

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:117125(ARRB1 antibody) at dilution of 1:100 (under 10x lens)

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IF, IP

## Dilutions


Recommended Dilution:

# beta Arrestin 1 antibody


Catalog Number: 108480



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:117125(ARRB1 antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using ARRB1 antibody Catalog No:117125 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-Arrestin (IP:Catalog No:117125, 4ug; Detection:Catalog No:117125 1:500) with mouse lung tissue lysate 4000ug.