

# Anti-GTPBP4/NOG1 Rabbit antibody



Catalog Number: 167880

## Product name

Rabbit IgG

Anti-GTPBP4/NOG1 Rabbit antibody

## Applications

## Specificity

WB, IHC-P, FC, IP

Human

## Dilutions

## Antibody description

WB: 1/1000

GTPBP4/NOG1 Rabbit monoclonal antibody

IHC: 1/50

## Preparation

FC: 1/20


Antigen: A synthetic peptide of human GTPBP4/NOG1

IP: 1/20

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Validations



## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

Western blot detection of GTPBP4/NOG1 in A549,HL-60 using GTPBP4/NOG1 Rabbit mAb(1:1000 diluted)

## Ig Type