

GSTA2 antibody

Catalog Number: 111180

Product name

GSTA2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GSTA2 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 26kd

Preparation

This antibody was obtained by immunization of GSTA2 recombinant protein (Accession Number: XM_011514532). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:50-1:200

IHC: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111180(GSTA2 antibody) at dilution of 1:100



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111180(GSTA2 antibody) at dilution of 1:100 (under 10x lens)