

Anti-Ezrin (2C12) Mouse antibody



Catalog Number: 168180

Product name

Anti-Ezrin (2C12) Mouse antibody

Specificity

Human, Mouse, Rat

Antibody description

Ezrin (2C12) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Ezrin protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Ezrin in MCF7,3T3,A549 and Rat Brain lysates using Ezrin mouse mAb(dilution 1:1000).Predicted band size:69kDa.Observed band size:81kDa.