

MECR antibody



Catalog Number: 112580

Product name

MECR antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MECR Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue. Positive WB detected in mouse heart tissue, human stomach tissue. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of MECR recombinant protein (Accession Number: NM_016011). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112580(MECR antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human heart using Catalog No:112580(MECR antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:112580(MECR antibody) at dilution of 1:100 (under 40x lens)