

USF1 antibody



Catalog Number: 116580

Product name

USF1 antibody

Immunogen

[Human USF1 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

USF1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, K-562 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of USF1 recombinant protein (Accession Number: NM_007122). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116580(USF1 Antibody) at dilution of 1:600



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:116580(USF1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)