

BRD2 antibody

Catalog Number: 108580

Product name

BRD2 antibody

IP: 1:200-1:1000

IF: 1:50-1:500


Specificity

Human; other species not tested.

Validations

Antibody description

BRD2 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 100-110 kDa




Preparation

This antibody was obtained by immunization of BRD2 recombinant protein (Accession Number: NM_001199455). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117227(BRD2 Antibody) at dilution of 1:600

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

IP Result of anti-BRD2 (IP:Catalog No:117227, 4ug; Detection:Catalog No:117227 1:300) with HeLa cells lysate 1200ug.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IP, IF

Dilutions

Recommended Dilution:

WB: 1:200-1:2000



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog

BRD2 antibody



Catalog Number: 108580

No:117227(BRD2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)