

Human CAMK1D/CAMKID (GST tag) recombinant protein



Catalog Number: 519580

General Information

Gene Name Synonym

Calcium/calmodulin-dependent protein kinase type 1D, CaM kinase I delta, CaMKI-like protein kinase, CAMK1D.CAMK1D, CaMK1 delta,

Protein Construction

A DNA sequence encoding the Human CAMK1D/CAMKID Met1-Lys385 is expressed with a GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB□ 150mM NaCl,pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human CAMK1D/CAMKID (GST tag) recombinant protein