

GIPC2 antibody



Catalog Number: 110980

Product name

GIPC2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GIPC2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 34 kDa

Preparation

This antibody was obtained by immunization of GIPC2 recombinant protein (Accession Number: NM_017655). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110980(GIPC2 antibody) at dilution of 1:300