

GIPC2 antibody



Catalog Number: 110980

Product name

Polyclonal

GIPC2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

GIPC2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 34 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:1000-1:10000

This antibody was obtained by immunization of GIPC2 recombinant protein (Accession Number: NM_017655). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110980(GIPC2 antibody) at dilution of 1:300

