

# MYO18A antibody



Catalog Number: 112980

## Product name

MYO18A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MYO18A Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, K-562 cells, Neuro-2a cells, NIH/3T3 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 230kd, 190kd

## Preparation

This antibody was obtained by immunization of MYO18A recombinant protein (Accession Number: NM\_001346767). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:1000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112980(MYO18A antibody) at dilution of 1:6000



IP Result of anti-MYO18A (IP:Catalog No:112980, 3ug; Detection:Catalog No:112980 1:300) with HeLa cells lysate 3000ug.