

# SLC16A4 Polyclonal Antibody



Catalog Number: 163980

## Product name

SLC16A4 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to SLC16A4

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 250-350 of human SLC16A4 (NP\_004687.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:1000 - 1:2000

## Validations



Western blot - SLC16A4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SLC16A4 Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.