

PNO1 antibody

Catalog Number: 113980

Product name

PNO1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PNO1 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, mouse spleen tissue. Positive IP detected in L02 cells. Observed molecular weight by Western-blot: 35kd

Preparation

This antibody was obtained by immunization of PNO1 recombinant protein (Accession Number: NM_020143). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:113980(PNO1 antibody) at dilution of 1:300



IP Result of anti-PNO1 (IP:Catalog No:113980, 3ug; Detection:Catalog No:113980 1:300) with L02 cells lysate 1500ug.