

# SBNO1 antibody

Catalog Number: 114980

## Product name

SBNO1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SBNO1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, mouse brain tissue, Raji cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 145-154 kDa

## Preparation

This antibody was obtained by immunization of SBNO1 recombinant protein (Accession Number: NM\_001167856). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114980(SBNO1 Antibody) at dilution of 1:600



IP Result of anti-SBNO1 (IP:Catalog No:114980, 3ug; Detection:Catalog No:114980 1:300) with HeLa cells lysate 1800ug.