

# PRSS2 antibody

Catalog Number: 114280

---

## Product name

PRSS2 antibody

## Specificity

Human,Mouse,Rat; other species not tested.

## Antibody description

PRSS2 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, HepG2 cells, mouse thymus tissue. Observed molecular weight by Western-blot: 26kd,35kd

## Preparation

This antibody was obtained by immunization of PRSS2 recombinant protein (Accession Number: BC030260). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:114280(PRSS2 antibody) at dilution of 1:400