

CHRNA5 antibody



Catalog Number: 109280

Product name

CHRNA5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CHRNA5 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, A431 cells, A549 cells, mouse bladder tissue, mouse stomach tissue. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of CHRNA5 recombinant protein (Accession Number: NM_000745). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109280(CHRNA5 antibody) at dilution of 1:500