

Anti-ING5 antibody



Catalog Number: 100680

Product name

Anti-ING5 antibody

Immunogen

[Human ING5 \(GST Tag\) recombinant protein](#)

Specificity

Human ING5

Antibody description

Rabbit polyclonal to ING5

Preparation

Produced in rabbits immunized with purified, recombinant Human ING5 (rh ING5; Q8WYH8-2; Glu36-Gln84). ING5 specific IgG was purified by Human ING5 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB


Dilutions

WB: 2-10 µg/mL

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human ING5. The detection limit for Human ING5 is approximately 0.00245 ng/well.

Validations



Items	Lanes	A
Sample (whole cell lysate)		HepG2
Sample Volume (µg/lane)		30
Gel		13% SDS PAGE reducing gel
Recommended Concentration		2-10 µg/ml
Secondary Antibody		Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.

ING5 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot