

# UPF1 Polyclonal Antibody



Catalog Number: 161680

## Product name

UPF1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to UPF1

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 270-370 of human UPF1 (NP\_002902.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations




Immunohistochemistry - UPF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using UPF1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - UPF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using UPF1 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - UPF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using UPF1 antibody at dilution of

# UPF1 Polyclonal Antibody



Catalog Number: 161680

---

1:100 (40x lens).