

H8N4 HA (His Tag) recombinant protein



Catalog Number: 503680

General Information

Protein Construction

A DNA sequence encoding the N-terminal segment (Met 1-Arg 344) of the influenza A hemagglutinin (A/pintail duck/Alberta/114/1979 (H8N4)) (ABB87729.1), termed as HA1, was expressed, fused with a C-terminal polyhistidine tag.

Organism

H8N4

Expression Host

Human Cells

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 18

Molecular Mass

The secreted recombinant influenza A H8N4 HA1 subunit (A/pintail duck/Alberta/114/1979 (H8N4)) comprises 338 amino acids and has a predicted

molecular mass of 38.6 kDa. As a result of glycosylation, it migrates as an approximately 50-60 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H8N4 (A/pintail duck/Alberta/114/1979) Hemagglutinin Protein (HA1 Subunit) (His Tag) SDS-PAGE