

PDHX antibody

Catalog Number: 113680

Product name

PDHX antibody

IP: 1:2000-1:20000

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description


PDHX Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, HEK-293 cells, SKOV-3 cells. Positive IP detected in HEK-293 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 54 kDa



Preparation

This antibody was obtained by immunization of PDHX recombinant protein (Accession Number: NM_003477). Purification method: Antigen affinity purified.

BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113680(PDHX antibody) at dilution of 1:5000



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113680(PDHX antibody) at dilution of 1:100 (under 40x lens)

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC, IP, WB

Dilutions

Recommended Dilution:

WB: 1:2000-1:20000



IP Result of anti-PDHX (IP:Catalog No:113680,

PDHX antibody



Catalog Number: 113680

4ug; Detection:Catalog No:113680 1:5000) with HEK-293 cells lysate 2000ug.