

Anti-ZNF217 Rabbit antibody



Catalog Number: 167680

Product name

Anti-ZNF217 Rabbit antibody

Specificity

Human

Antibody description

ZNF217 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human ZNF217

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/2000-1/10000

Validations



Western blot detection of ZNF217 in A549,HL-60,U251 cell lysates using ZNF217 Rabbit mAb(1:1000 diluted).Predicted band size:115kDa.Observed band size:115kDa.