

# Mouse B7-H5 / GI24 / VISTA (Fc Tag) recombinant protein



Catalog Number: 500380

## General Information

### Protein Construction

A DNA sequence encoding the mouse C10orf54 (NP\_083008.1) (Met1-Ala191) was expressed with the Fc region of human IgG1 at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Phe 33

### Molecular Mass

The recombinant mouse C10orf54 consists 397

amino acids and predicts a molecular mass of 44.5 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse B7-H5 / GI24 / VISTA Protein (Fc Tag)