

Mouse BMP-4/BMP-2B (Fc Tag) recombinant protein



Catalog Number: 502380

General Information

Gene Name Synonym

Bone morphogenetic protein 2B

Protein Construction

A DNA sequence encoding the mouse BMP4 (NP_031580.2) (Ser293-Arg408) was expressed with the Fc region of rabbit IgG at the N-terminus..

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ser

Molecular Mass

The recombinant mouse BMP4 consists 365 amino acids and predicts a molecular mass of 40.6 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse BMP4 / BMP-4 Protein (Fc Tag)