

MAGEA4 antibody

Catalog Number: 112380

Product name

MAGEA4 antibody

Specificity

Human; other species not tested.

Antibody description

MAGEA4 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, A375 cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 35-40 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001011548). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:112380(MAGEA4 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:112380(MAGEA4 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:112380(MAGEA4 antibody) at dilution of 1:50 (under 40x lens)