

# IL12B antibody



Catalog Number: 107380

## Product name

IL12B antibody

## Specificity

Human, Pig; other species not tested.

## Antibody description

IL12B Mouse Monoclonal antibody. Positive WB detected in human liver tissue. Observed molecular weight by Western-blot: 37 kDa

## Preparation

This antibody was obtained by immunization of IL12B recombinant protein (Accession Number: NM\_002187). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107380(IL12B Antibody) at dilution of 1:1000