

Anti-Bax (2B12) Mouse antibody



Catalog Number: 168380

Product name

Anti-Bax (2B12) Mouse antibody

Ig Type

Mouse IgG2a

Specificity

Human, Monkey

Applications

WB

Antibody description

Bax (2B12) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Purified recombinant human Bax protein fragments expressed in E.coli

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of Bax in COS7,Raji,CEM and HeLa cell lysates using Bax mouse mAb (1:1000 diluted). Predicted band size:20KDa. Observed band size:20KDa.