

Collagen Type IV antibody



Catalog Number: 109380

Product name

Collagen Type IV antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Collagen Type IV Rabbit Polyclonal antibody.
Positive WB detected in mouse heart tissue,
human heart tissue, mouse kidney tissue.
Observed molecular weight by Western-blot: 50
kDa

Preparation

This antibody was obtained by immunization of
Collagen Type IV recombinant protein (Accession
Number: NM_001287760). Purification method:
Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse heart tissue were subjected to SDS PAGE
followed by western blot with Catalog
No:109380(COL4A6 antibody) at dilution of
1:1000