

# Anti-KDR antibody



Catalog Number: 142780

## Product name

Anti-KDR antibody

## Immunogen

[Mouse VEGFR2/Flk-1/CD309/KDR \(His Tag\) recombinant protein](#)

## Specificity

Mouse VEGFR2 / Flk-1 / CD309 / KDR

## Antibody description

Rabbit Monoclonal to Mouse KDR

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse VEGFR2 / Flk-1 / CD309 / KDR and conjugated with APC under optimum conditions, the unreacted APC was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from

prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

FCM

## Dilutions

FCM

## Validations



VEGFR2 / Flk-1 / CD309 / KDR Antibody (APC), Rabbit MAb, Flow Cytometry