

SYNCRIP antibody

Catalog Number: 115780

Product name

SYNCRIP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SYNCRIP Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HEK-293 cells, Jurkat cells. Positive IP detected in A431 cells. Observed molecular weight by Western-blot: 62-74 kDa

Preparation

This antibody was obtained by immunization of SYNCRIP recombinant protein (Accession Number: NM_001159674). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:115780(SYNCRIP antibody) at dilution of 1:500



IP Result of anti-SYNCRIP (IP:Catalog No:115780, 4ug; Detection:Catalog No:115780 1:500) with A431 cells lysate 1200ug.