

VPS33A antibody

Catalog Number: 116780

Product name

VPS33A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VPS33A Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in human brain tissue, HepG2 cells, human cerebellum tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of VPS33A recombinant protein (Accession Number: NM_022916). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116780(VPS33A antibody) at dilution of 1:500



IP Result of anti-VPS33A (IP:Catalog No:116780, 4ug; Detection:Catalog No:116780 1:600) with HepG2 cells lysate 1600ug.