

DDX10 antibody

Catalog Number: 109780

Product name

DDX10 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DDX10 Rabbit Polyclonal antibody. Positive IHC detected in mouse testis tissue, human kidney tissue. Positive WB detected in SMMC-7721 cells, HeLa cells. Observed molecular weight by Western-blot: 101 kDa

Preparation

This antibody was obtained by immunization of DDX10 recombinant protein (Accession Number: NM_004398). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with Catalog No:109780(DDX10 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:109780(DDX10 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded mouse testis tissue slide using Catalog No:109780(DDX10 Antibody) at dilution of 1:50

DDX10 antibody



Catalog Number: 109780

(under 40x lens)