

# DNAJB12 antibody

Catalog Number: 110090

## Product name

DNAJB12 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DNAJB12 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse uterus tissue, human brain tissue, human kidney tissue, human liver tissue. Observed molecular weight by Western-blot: 42 kDa

## Preparation

This antibody was obtained by immunization of DNAJB12 recombinant protein (Accession Number: BC011812). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse uterus tissue were subjected to SDS PAGE followed by western blot with Catalog No:110090(DNAJB12 antibody) at dilution of



IP Result of anti-DNAJB12 (IP:Catalog No:110090, 4ug; Detection:Catalog No:110090 1:500) with mouse brain tissue lysate 4000ug.