

CDAN1 antibody

Catalog Number: 109090

Product name

CDAN1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CDAN1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, K-562 cells. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 130 kDa

Preparation

This antibody was obtained by immunization of CDAN1 recombinant protein (Accession Number: NM_138477). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109090(CDAN1 antibody) at dilution of 1:500



IP Result of anti-CDAN1 (IP:Catalog No:109090, 4ug; Detection:Catalog No:109090 1:500) with K-562 cells lysate 3200ug.