

# Anti-KIFAP3 Rabbit antibody



Catalog Number: 167890

## Product name

Anti-KIFAP3 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

KIFAP3 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human KIFAP3

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/1000

IHC: 1/50

ICC/IF: 1/20

## Validations



Western blot detection of KIFAP3 in A549,HL-60,C6 cell lysates using KIFAP3 Rabbit mAb(1:1000 diluted).Predicted band size:91kDa.Observed band size:91kDa.