

# ADRA1B-Specific antibody



Catalog Number: 107890

## Product name

ADRA1B-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ADRA1B-Specific Rabbit Polyclonal antibody.  
Positive IHC detected in mouse brain tissue,  
human heart tissue.

## Preparation

This antibody was obtained by immunization of  
Peptide (Accession Number: NM\_000679).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH  
7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:50-1:500

## Validations



Immunohistochemistry of paraffin-embedded  
mouse brain tissue slide using Catalog  
No:107890(ADRA1B-Specific Antibody) at dilution  
of 1:100 (under 10x lens). heat mediated antigen  
retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded  
mouse brain tissue slide using Catalog  
No:107890(ADRA1B-Specific Antibody) at dilution  
of 1:100 (under 40x lens). heat mediated antigen  
retrieved with Tris-EDTA buffer(pH9).