

Human PIN4 (His tag) recombinant protein



Catalog Number: 519890

General Information

Gene Name Synonym

Peptidyl-prolyl cis-trans isomerase NIMA-interacting 4, Parvulin-14, Parvulin-17, Peptidyl-prolyl cis-trans isomerase Pin4, Peptidyl-prolyl cis/trans isomerase EPVH, Rotamase Pin4, PIN4,

Protein Construction

A DNA sequence encoding the Human PIN4 Met1-Lys156 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PIN4 (His tag) recombinant protein