

# LDOC1 antibody



Catalog Number: 112190

## Product name

LDOC1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

LDOC1 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 23 kDa

## Preparation

This antibody was obtained by immunization of LD<sup>o</sup>C1 recombinant protein (Accession Number: NM\_012317). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112190(LDoc1 Antibody) at dilution of 1:350