

# NIPSNAP3B antibody



Catalog Number: 113190

## Product name

NIPSNAP3B antibody

## Specificity

Human; other species not tested.

## Antibody description

NIPSNAP3B Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IHC detected in human lymphoma tissue. Observed molecular weight by Western-blot: 28kd

## Preparation

This antibody was obtained by immunization of NIPSNAP3B recombinant protein (Accession Number: NM\_018376). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113190(NIPSNAP3B antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:113190(NIPSNAP3B antibody) at dilution of 1:50 (under 10x lens)