

ACADSB Polyclonal Antibody



Catalog Number: 164190

Product name

ACADSB Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ACADSB

Preparation

Antigen: A synthetic peptide of human ACADSB

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - ACADSB Polyclonal Antibody

Western blot analysis of extracts of 293 cells, using ACADSB antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.