

Cyclin D1 antibody



Catalog Number: 107190

Product name

Cyclin D1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Cyclin D1 Mouse Monoclonal antibody. Positive WB detected in NIH/3T3 cells, A549 cells, HeLa cells, SH-SY5Y cells. Observed molecular weight by Western-blot: 36 kDa

Preparation

This antibody was obtained by immunization of Cyclin D1 recombinant protein (Accession Number: NM_053056). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1000-10000

Validations



Western blot of Cyclin D1 in HepG2, SW 1990 and NIH/3T3 cell lines with Catalog No:107190 at dilution of 1:10000