

Anti-GSK3 beta (9F9) Mouse antibody



Catalog Number: 168190

Product name

Anti-GSK3 beta (9F9) Mouse antibody

Specificity

Human, Rat, Mouse

Antibody description

GSK3 beta (9F9) Mouse monoclonal antibody

Preparation

Antigen: Synthetic peptide conjugated to KLH.

Formulation

PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, IHC-P

Dilutions


WB: 1/1;000-2;000

IHC: 1/100-200


Validations



Immunohistochemical analysis of paraffin-embedded Human Stomach Carcinoma Tissue using GSK3 β Mouse mAb diluted at 1:200.



Western blot detection of GSK3 β in Goat muscle, CHO-K1, C6, HeLa, 3T3 cell lysates using GSK3 β (9F9) Mouse mAb(1:1000 diluted). Predicted band size:46KDa. Observed band size:46KDa.




Western blot analysis of 1) HeLa Cell Lysate, 2) 3T3 Cell Lysate, 3) Rat Brain Tissue Lysate using GSK3 β Mouse mAb diluted at 1:1000.

Anti-GSK3 beta (9F9) Mouse antibody



Catalog Number: 168190



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using GSK3 β Mouse mAb diluted at 1:200.