

# ARHGEF9-Specific antibody



Catalog Number: 108190

## Product name

ARHGEF9-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ARHGEF9-Specific Rabbit Polyclonal antibody.  
Positive WB detected in human brain tissue, HeLa cells. Observed molecular weight by Western-blot: 56 kDa, 61 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_015185).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108190(ARHGEF9-Specific antibody) at



dilution of 1:300

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108190(ARHGEF9-Specific antibody) at dilution of 1:200