

CERCAM antibody



Catalog Number: 109190

Product name

Polyclonal

CERCAM antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

CERCAM Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 68kd,48kd

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of CERCAM recombinant protein (Accession Number: BC011811). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109190(CERCAM antibody) at dilution of 1:1500