

# CYBA Polyclonal Antibody



Catalog Number: 160590

---

## Product name

CYBA Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to CYBA

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 126-195 of human CYBA (NP\_000092.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:100

## Validations



Immunofluorescence - CYBA Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using CYBA antibody at dilution of 1:100 (40x lens).  
Blue: DAPI for nuclear staining.