

BRIT1 antibody



Catalog Number: 108590

Product name

BRIT1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

BRIT1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, L02 cells, MDA-MB-453s cells, mouse ovary tissue. Positive IHC detected in human prostate cancer tissue, human kidney tissue. Observed molecular weight by Western-blot: 66 kDa, 100 kDa

Preparation

This antibody was obtained by immunization of BRIT1 recombinant protein (Accession Number: NM_001172574). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117237(BRIT1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:117237(BRIT1 antibody) at dilution of 1:50 (under 10x lens)