

IGHMBP2 Polyclonal Antibody



Catalog Number: 162990

Product name

IGHMBP2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to IGHMBP2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 540-720 of human IGHMBP2 (NP_002171.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations



Immunohistochemistry - IGHMBP2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using IGHMBP2 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - IGHMBP2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using IGHMBP2 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - IGHMBP2 Polyclonal Antibody

IGHMBP2 Polyclonal Antibody



Catalog Number: 162990

Immunohistochemistry of paraffin-embedded human mammary cancer using IGHMBP2 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - IGHMBP2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver using IGHMBP2 antibody at dilution of 1:100 (40x lens).