

# Podoplanin antibody



Catalog Number: 113990

## Product name

Podoplanin antibody

## Specificity

Human; other species not tested.

## Antibody description

Podoplanin Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, human placenta tissue. Observed molecular weight by Western-blot: 43 kDa

## Preparation

This antibody was obtained by immunization of Podoplanin recombinant protein (Accession Number: BC022812). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113990(PDPN antibody) at dilution of 1:1000